DK

Axis® Homocysteine EIA

REF FHCY100

AXIS-SHIELD

Axis-Shield Diagnostics Ltd.
The Technology Park
Dundee DD2 1XA
United Kingdom
Tel: +44 (0) 1382 422000
Fax: +44 (0) 1382 422088
www.axis-shield.com

REF

Produktnummer

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IVD

Til *in vitro*-diagnostisk brug

Udløpsdato

LOT

Partnummer

Opbevares ved 2-8 °C

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
Axis® Homocysteine EIA - University of Bergen HPLC

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Referencer

1. Malinow, MR. Plasma homocyst(e)ine and arterial occlusive diseases: A mini-review. *Clin Chem* 1995;41:173-176.
2. Ueland PM. Homocysteine species as components of plasma redox thiol status. *Clin Chem* 1995;41:340-342.
3. Perry IJ, Refsum H, Morris RW, Ebrahim SB, Ueland PM, Shaper AG. Prospective study of serum total homocysteine concentration and risk of stroke in middle-aged British men. *The Lancet* 1995;346:1395-1398.
4. Finkelstein JD. Methionine metabolism in mammals. *J Nutr Biochem* 1990;1:228-237.
5. Mudd SH, Levy HL, Skovby F. Disorder of transsulfuration. In: C.R. Scriver, A.L. Beaudet, W.S. Sly, et al. eds. *The metabolic basis of inherited disease* New York: McGraw-Hill; 1995:1297-1327.
6. Clarke R, Daly L, Robinson K, Naughten E, Cahalane S, Fowler B, Graham I. Hyperhomocysteinemia: An independent risk factor for vascular disease. *N Engl J Med* 1991;324:1149-1155.
7. Deloughery TG, Evans A, Sadeghi A, McWilliams J, Henner WD, Taylor LM Jr, Press RD. Common mutation in methylenetetrahydrofolate reductase. Correlation with homocysteine metabolism and late-onset vascular disease. *Circulation* 1996;94:3074-3078.
8. Schmitz C, Lindpaintner K, Verhoef P, Gaziano JM, Buring J. Genetic polymorphism of methylenetetrahydrofolate reductase and myocardial infarction. *Circulation* 1996;94:1812-1814.
9. Boushey CJ, Beresford SAA, Omenn GS, Motulsky AG. A Quantitative Assessment of Plasma Homocysteine as a Risk Factor for Vascular Disease. *JAMA* 1995;274:1049-1057.
10. Allen RH, Stabler SP, Savage DG, Lindenbaum J. Diagnosis of cobalamin deficiency I: Usefulness of serum methylmalonic acid and total homocysteine concentrations. *Am J Hematol* 1990;34:90-98
11. Stabler SP, Lindenbaum J, Allen RH. The Use of homocysteine and other metabolites in the specific diagnosis of vitamin B-12 deficiency. *J Nutr* 1996;126: 12665-12725.
12. Guttormsen AB, Svarstad E, Ueland PM, Refsum H. Elimination of homocysteine from plasma in subjects with end-stage renal failure. *Irish J Med Sci* 1995;164 (Suppl. 15):8-9
13. Bostom AG, Lathrop L. Hyperhomocysteinemia in end-stage renal disease: Prevalence, etiology, and potential relationship to arteriosclerotic outcomes. *Kidney Int* 1997;52:10-20.
14. Desouza C, Keebler M, McNamara DB, Fonseca V. Drugs affecting homocysteine metabolism. *Review. Drugs* 2002;62 (4):605-616.
15. Refsum H, Ueland PM. Clinical significance of pharmacological modulation of homocysteine metabolism. *Review. Trends Pharmacol Sci* 1990;11:411-416
16. Frantzen F, Faaren AL, Alfheim I, Nordhei AK. An enzyme conversion immunoassay for determining total homocysteine in plasma or serum. *Clin Chem* 1998;44:311-316.
17. Sundrehagen E. Axis Biochemicals ASA. Enzymatic assay for homocysteine and a kit therefor. EP 623174/US5631127.
18. US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. HHS Publication No. (CDC) 93-8395. Washington, DC: US Government Printing Office, May 1999 (4th edition).
19. Ubbink JB, Vermaak Hayward WJ, van der Merwe A, Becker PJ. The effect of blood sample aging and food consumption on plasma total homocysteine levels. *Clinica Chimica Acta* 1992;207:119-128.
20. Guttormsen AB, Schneede J, Fiskerstrand T, Ueland PM, Refsum H. Plasma concentrations of homocysteine and other aminothioli compounds are related to food intake in healthy human subjects. *J Nutr* 1994;124:1934-1941.
21. Ueland PM, Refsum H, Stabler SP, Malinow R, Andersson A, Allen RH. Total homocysteine in plasma or serum: Methods and clinical applications. *Clin Chem* 1993;39:1764-1779.
22. Nehler MR, Taylor LM Jr, Porter JM. Homocysteine as a risk factor for atherosclerosis: A review. *Cardiovasc Pathol* 1997;6:1-9.
23. Lussier-Cacan S, Xhignesse M, Piolot A, Selhub J, Davignon J, Genest J Jr. Plasma total homocysteine in healthy subjects: sex-specific relation with biological traits. *Am J Clin Nutr* 1996; 64:587-593.
24. Clarke R, Woodhouse P, Ulvik A, Frost C, Sherliker P, Refsum H, Ueland PM, Khaw K- T. Variability and determinants of total homocysteine concentrations in plasma in an elderly population. *Clin Chem* 1998;44:102-107.
25. Jacques PF, Selhub J, Bostom AG, Wilson PWF, Rosenberg IH. The effect of folic acid fortification on plasma folate and total homocysteine concentrations. *N Engl J Med* 1999;340:1449-1454.
26. Lawrence JM, Pettitt DB, Watkins M, Umekubo MA. Trends in serum folate after food fortification. *Lancet* 1999;354:915-916.
27. Fiskerstrand T, Refsum H, Kvalheim G, Ueland PM, Homocysteine and other thiols in plasma and urine: Automated determination and sample stability. *Clin Chem* 1993;39:263-271.